

**TIM TRIO**

**Basic Scan Procedure**

# Scan Procedure

- The scan procedure has not changed much with the TIM TRIO system. You may notice that there are some pop-menus that you have not seen on the previous TRIO system.
- Also some of the pop-up menus, although are the same as the TRIO, it just looks different. The following pages goes through the basic scan procedure.
- Always remember to register your subject first before you position and landmark. If you don't the scan location on the LCD display will change from your original landmark without notice.

- Register... Num 0
- Browser... Num .
- Schedule...
- Close Patient
- Save As...
- Copy Selection
- Reset Table Position
- Voice Output Properties...
- Filming Layout...
- Film Task Status...
- Film Preview...
- Expose Film Task Ctrl+P
- Copy to Film Sheet
- Search
- Search Selected...

Click patient and select register.

Exam

Viewing

Filming

3D

<Please register new patient>




Open Apply Cancel

Region:

Exams:

Programs:



Stop Continue Skip

Program



Complete all bold fields.

Click on "Exam".

**Patient Registration**

<b>PATIENT</b>	<b>Last name</b>	test
	<b>First name</b>	
	<b>Title</b>	
	<b>Patient ID</b>	111
	<b>Date of Birth</b>	4/4/1987 [M/d/yyyy]
	<b>Sex</b>	<input checked="" type="radio"/> Male <input type="radio"/> Female <input type="radio"/> Other
	<b>Age</b>	22 Years
	<b>Height</b>	[ft'in"]
	<b>Weight</b>	180.00 [lb]
	<b>Additional info</b>	
<input type="button" value="Details..."/>		
<b>HOSPITAL</b>	<b>Referring physician</b>	
	<b>Requesting physician</b>	
	<b>Admission ID</b>	
<input type="button" value="Preregister"/> <input type="button" value="Exam"/> <input type="button" value="Search"/> <input type="button" value="Cancel"/> <input type="button" value="Help"/>		

<b>PROCEDURE</b>	<b>Accession No</b>		
	<b>Request ID</b>		
	<b>Requested procedure(s)</b>	Warmup warmup	
	<b>Study</b>	Warmup warmup	
	<b>Study comment</b>		
	<b>Patient position</b>	Head First - Supine	
	<input type="button" value="Details..."/>		
	<b>INSTITUTION</b>	<b>Institution name</b>	
		<b>1. Performing physician</b>	
		<b>1. Operator</b>	
<input type="button" value="Details..."/>			

Exam

Viewing

Filming

3D



Stop

Continue

Skip

Program

Click on "Confirm".

**Patient Registration**

<b>PATIENT</b>	Last name	test
	First name	
	Title	
	Patient ID	111
	Date of Birth	4/4/1987 [M/d/yyyy]
	Sex	<input checked="" type="radio"/> Male <input type="radio"/> Female <input type="radio"/> Other
	Age	
	Height	
	Weight	
	Additional info	

<b>PROCEDURE</b>	Accession No	
	Request ID	
	Requested procedure(s)	<input type="checkbox"/> Warmup warmup
	Study	Warmup warmup

**Registration Confirmation**

Patient name: test

Weight: [lb] 180.00

Patient Position: Head First - Supine

<b>HOSPITAL</b>	Referring physician	
	Requesting physician	
	Admission ID	

<b>STATION</b>	1. Performing physician	
	1. Operator	

Exam

Viewing

Filming

3D



Stop

Continue

Skip

Program

Exam

Viewing

Filming

3D

Select sequence and drag and drop to queue or click "<<" button

test 4/4/1987

Open Apply Cancel

Region: Warmup

Exams: warmup

Programs: Warm\_up001

	localizer	00:09
2	epi2d_bold	04:52
3	epi2d_bold	04:52

<<

Program

Σ 00:00 Σ 09:53

Stop Continue Skip

0%

Exam

Viewing

Filming

3D

Click on continue if the following pop-up appears.

**Warning**

 **Automatic Table Move**

Further measurements could automatically move the table.  
Ensure that the patient is prepared for table movement.

Table Move Notification

Continue Cancel Help

test

Open Apply

1 localizer

Exams: warmup

Programs:

Warm_up001	1	localizer	00:09
	2	 epi2d_bold	04:52
	3	 epi2d_bold	04:52

<<

Σ 00:00

Σ 09:53


Program

Preparing...

Stop Continue Skip

SAR=NM





Click on "Apply" to start scan. Images will appear in viewports.

test 4/4/1987

Scan Apply Cancel

1 localizer

Region: Warmup

Exams: warmup

Programs:			
Warm_up001	1	localizer	00:09
	2	epi2d_bold	04:52
	3	epi2d_bold	04:52


Σ 00:00

Stop Continue Skip

Program

Σ 09:53

5%



Select next sequence to be scanned and drop in queue. Click on "Open" to activate slice grid.

test 4/4/1987

Open Apply Cancel

1	localizer	
2	epi2d_bold	04:52

Σ 04:52

Stop Continue Skip

Region: Warmup

Exams: warmup

Programs: Warm\_up001

1	localizer	00:09
2	epi2d_bold	04:52
3	epi2d_bold	04:52




<<

Program

Σ 09:53

Exam  
Viewing  
Filming  
3D

TP 0  
SP F1.4  
FoV 280\*280  
Tra



Exam

Viewing

Filming

3D

Position slices as needed. Click on "Apply"

test 4/4/1987 TA: 4:52 PM: REF PAT: Off Voxel size: 3.4x3.4x4.0 mm Rel. SNR: 1.00 : epfid

Scan Apply Cancel

1 localizer

2 epi2d\_bold

Σ 00:00

Stop Continue Skip

Slice group 1

Slices 34

Dist. factor 0 %

Position R1.4 A31.2 F1.2

Orientation T > C-9.8

Phase enc. dir. A >> P

Phase oversampling 0 %

FoV read 220 mm

FoV phase 100.0 %

Slice thickness 4.0 mm

TR 2000 ms

TE 25 ms

Averages 1

Concatenations 1


Filter None

Coil elements PH

Program Routine Contrast Resolution Geometry System Physio BOLD Sequence

Parent Applications Transfer Edit Queue Protocol View Image Tools Evaluation Scroll System Add-On Options Help


test 111  
1 IMA 1 / 3  
9/17/2009



test 111  
1 IMA 2 / 3  
9/17/2009



test 111  
1 IMA 3 / 3  
9/17/2009



Exam

Viewing

Filming

3D

Exam pause pop-up appears. Click "Continue" to begin scan.

M/ND  
localizer

FoV 280\*280  
Sag

FoV 280\*280  
Cor

localizer

TP 0  
SP F1.4  
FoV 280\*280  
Tra

test 4/4/1987

Open Apply Cancel

1	localizer	
2	epi2d_bold	00:34

**Exam paused**

2 epi2d\_bold

Patient has contrast agent

Continue Close Skip Help

Region: Warmup		
Exams: warmup		
Programs: Warm_up001		
1	localizer	00:09
2	epi2d_bold	04:52
3	epi2d_bold	04:52


Σ 00:34

Σ 09:53

Stop Continue Skip

Waiting for user to continue. Stimu=NM SAR=NM

9/17/2009 8:54:24 AM



Exam

Viewing

Filming

3D

Select next sequence to scan and drop in queue. Click "Open" to activate slice grid.

TP 0  
SP F1.4  
FoV 280\*280  
Tra

Buttons: **Open** (highlighted with red arrow), **Apply**, **Cancel**

1	localizer	00:09
2	epi2d_bold	04:52
3	epi2d_bold	04:52

Σ 04:52


Buttons: **Stop**, **Continue**, **Skip**

Region: Warmup  
Exams: warmup  
Programs: Warm\_up001

1	localizer	00:09
2	epi2d_bold	04:52
3	epi2d_bold	04:52

Program

Σ 09:53



To copy: Highlight sequence to copy. In this case it is sequence # 2.




test 4/4/1987 TA: 4:52 PM: REF PAT: Off Voxel size: 3.4x3.4x4.0 mm Rel. SNR: 1.00 : epid

1	localizer	
2	<b>epi2d_bold</b>	
3	epi2d_bold	

Σ 00:00

Slice group	1	FoV read	220	mm
Slices	34	FoV phase	100.0	%
Dist. factor	0	Slice thickness	4.0	mm
Position	R1.4 A31.2 F1.2	TR	2000	ms
Orientation	T > C-9.8	TE	25	ms
Phase enc. dir.	A >> P	Averages	1	
Phase oversampling	0	Concatenations	1	
		Filter	None	
		Coil elements	PH	

Program Routine Contrast Resolution Geometry System Physio BOLD Sequence



Exam

Viewing

Filming

3D

Right click and select "Copy Parameter"

- Open
- View Protocol...
- Rerun
- Complete
- Copy Parameter...
- Update Copy Reference...
- Append
- Move Up
- Move Down
- Load Series
- Cut
- Copy
- Paste
- Change Procedure Step...
- Change Body Part Examined...
- Set Image Comment...
- Delete
- Clear All
- Properties...

M/ND  
localizer

test

Scan

1 localizer

2 epi2d bold

3 epi2d bold

TP 0  
SP A67.6  
FoV 280\*280  
Cor

M/ND  
localizer

TP 0  
SP F1.4  
FoV 280\*280  
Tra

TA: 4:52 PM: REF PAT: Off Voxel size: 3.4x3.4x4.0 mm Rel. SNR: 1.00 : epfid

Slice group	1	FoV read	220	mm
Slices	34	FoV phase	100.0	%
Dist. factor	0	Slice thickness	4.0	mm
Position	R1.4 A31.2 F1.2	TR	2000	ms
Orientation	T > C-9.8	TE	25	ms
Phase enc. dir.	A >> P	Averages	1	
Phase oversampling	0	Concatenations	1	
		Filter	None	
		Coil elements	PH	

Stop Continue Skip

Parent Applications Transfer Edit Queue Protocol View Image Tools Evaluation Scroll System Add-On Options Help

test 111  
1 IMA 1 / 3  
9/17/2009

test 111  
1 IMA 2 / 3  
9/17/2009

test 111  
1 IMA 3 / 3  
9/17/2009

Exam  
Viewing  
Filming  
3D

M/ND  
localizer

test  
Scan Apply

1 localizer  
2 epi2d\_bold  
3 epi2d\_bold

Orientation: T > C-9.8  
Phase enc. dir: A >> P  
Phase oversampling: 0 %

×3.4×4.0 mm Rel. SNR: 1.00 : epi2d

FoV read: 220 mm  
FoV phase: 100.0 %  
Slice thickness: 4.0 mm  
TR: 2000 ms  
TE: 25 ms  
Averages: 1  
Concatenations: 1  
Filter: None  
Coil elements: PH

TP 0  
SP F1.4  
FoV 280\*280  
Tra

Σ 00:00

Stop Continue Skip

Program Routine Contrast Resolution Geometry System Physio BOLD Sequence

Waiting for slice positioning. SAR=NM 9/17/2009 9:00:19 AM




**Copy Parameter Groups**

- Center of slice groups & sat. regions
- Slices & sat. regions
- Slices
- Sat. regions
- Measurement parameters
- Adjust volume
- Everything
- Table position
- Navigators

Copy phase encoding direction  
 Steps

OK Cancel Help

Select "Center of slice groups & sat regions"  
Click "OK".



Exam

Viewing

Filming

3D

Click "Apply" to start scan

test 4/4/1987 TA: 4:52 PM: FIX PAT: Off Voxel size: 3.4x3.4x4.0 mm Rel. SNR: 1.00 : epfid

Scan Apply Cancel

- 1 localizer
- 2 epi2d\_bold
- 3 epi2d\_bold

Σ 00:00

Stop Continue Skip

Slice group 1 Slices 34 Dist. factor 0 Position R1.4 A31.2 F1.2 Orientation T > C-9.8 Phase enc. dir. A >> P Phase oversampling 0

FoV read 220 mm FoV phase 100.0 % Slice thickness 4.0 mm TR 2000 ms TE 25 ms Averages 1 Concatenations 1 Filter None Coil elements PH


Program Routine Contrast Resolution Geometry System Physio BOLD Sequence

Waiting for slice positioning.

SAR=NM



9/17/2009 9:01:35 AM



Exam pause pop-up appears. Click "Continue" to begin scan. For subsequent scans follow the same procedure.

1	localizer	00:09
2	epi2d_bold	04:52
3	epi2d_bold	00:34

**Exam paused**

3 epi2d\_bold

Patient has contrast agent

Continue Close Skip Help

Region: Warmup

Exams: warmup

Programs:

Warm_up001	1	localizer	00:09
	2	epi2d_bold	04:52
	3	epi2d_bold	04:52

Program

Patient Applications Transfer Edit Queue Protocol View Image Tools Evaluation Scroll System Add-On Options Help

Register... Num 0  
 Browser... Num .  
 Schedule...  
**Close Patient**  
 Save As...  
 Copy Selection  
 Reset Table Position  
 Voice Output Properties...  
 Filming Layout...  
 Film Task Status...  
 Film Preview...  
 Expose Film Task Ctrl+P  
 Copy to Film Sheet  
 Search  
 Search Selected...

**Exam**

**Viewing**

**Filming**

**3D**

When study is complete, Click "Patient" and select "Close Patient" to end study.

M/ND  
localizer

test 4/4/1987

1 localizer  
2 epi2d\_bold  
3 epi2d\_bold

Σ 00:00

Stop Continue Skip

TP 0  
SP F1.4  
FoV 280\*280  
Tra

Region: Warmup

Exams: warmup

Programs: Warm\_up001

1	localizer		00:09
2	epi2d_bold		04:52
3	epi2d_bold		04:52

Σ 09:53

Program

Waiting for scan instructions.

SAR=NM

9/17/2009 9:04:55 AM