

BD Alaris IV Infusion System Skills Checklist #DAHS-NSCBD18-ALARIS

Page 1 of 4

| | | |
|------------------|---|-----------------------|
| Name: | | Employee ID #: |
| Unit: | | Title: |
| Due Date: | New hire: prior to end of unit orientation period: ____/____/____. | |
| | Current Staff: | |

These skills will be considered complete when all below performance criteria are completed and pages 1, 2, 3 and 4 have been scanned and emailed to: hs-cppn@ucdavis.edu

Not all skills are applicable to all Nursing areas – if not applicable mark as N/A

BD Alaris IV Infusion System Skills Checklist #DAHS-NSCBD18-ALARIS

| References: | Date Completed (or N/A) | Verifier Initials |
|--|--------------------------------|--------------------------|
| 1. UC Davis Health Policy 13056: Parenteral Infusion Pump Use | | |
| 2. UC Davis Health Policy 13016: Intravenous Patient Controlled Analgesia Infusion - Adult/Pediatric | | |
| 3. UC Davis Health Policy 13033: Administration of Adult and Pediatric IV Medications | | |

Complete the assigned Alaris eLearning modules in UC Learning

BD Alaris IV Infusion System policies and procedures reviewed

Alaris™ Pump module

Demonstrate Pump Setup

- The patient's heart level should be in line with [CHANNEL SELECT] key
- Closes the administration set roller clamp when the safety clamp is open, to prevent free flow
- Does not use needles or blunt cannulas to access a SmartSite™ Needle-Free Valve
- Swabs the SmartSite™ Needle-Free Valve with a sterile 70% isopropyl alcohol wipe prior to any connection

Demonstrate System Start Up and Operation

- Understanding of what happens when [NEW PATIENT] is selected
- Understanding of the Patient Care Profile and how to change it

Demonstrate Programming with Guardrails™ Safety Software

- Programming a primary infusion on the Alaris™ Pump module
- Responding to a Guardrails™ Soft or Hard Limit alarm with audio alerts and visual prompts
- Programming an intermittent infusion on the Alaris™ Pump module
- Programming a Volume/Duration infusion on the Alaris™ Pump module
- Use of the "RESTORE" feature (previous programming, VTBI, bolus)
- Programming a medication bolus and describing the "Rapid Bolus" infusion feature
- Pausing an infusion by pressing the [PAUSE] hard key on the pump module and the PC unit
- The appropriate head height differential when hanging a 2° medication bag, or a 2° medication bottle

Demonstrate Basic Programming Without Guardrails™ Safety Software

- Programming of a Basic Infusion. Verbalize safety concerns when this mode is used

Alaris™ Syringe module

Demonstrate Syringe Module Setup

The patient's heart level should be in line with [CHANNEL SELECT] key.

- Priming the set using the Syringe Channel Option feature "Prime Set with Syringe." (*Infant and Child Only*)
- Proper priming technique when using an administration set with Pressure Sensing Disc tubing. (*Infant and Child Only*)

BD Alaris IV Infusion System Skills Checklist #DAHS-NSCBD18-ALARIS

Page 2 of 4

| | | |
|------------------|---|-----------------------|
| Name: | | Employee ID #: |
| Unit: | | Title: |
| Due Date: | New hire: prior to end of unit orientation period: ____ / ____ / ____. | |
| | Current Staff: | |

These skills will be considered complete when all below performance criteria are completed and pages 1, 2, 3 and 4 have been scanned and emailed to: hs-cppn@ucdavis.edu

| Alaris™ Syringe module (Continued) | Date Completed (or N/A) | Verifier Initials |
|------------------------------------|-------------------------|-------------------|
|------------------------------------|-------------------------|-------------------|

- Clamping the tubing after priming to prevent uncontrolled flow
- Loading and unloading a syringe into the Alaris Syringe module
- Correct selection of syringe manufacturer and size

- Recommend measures to help reduce start-up delays. (*Infant and Child Only*)

Demonstrate Basic Programming **Without** Guardrails™ Safety Software

- Programming of a Basic Infusion. Verbalize safety concerns when this mode is used

Demonstrate and Verbalize Measures to help Reduce Start-Up Delays (*Infant and Child Only*)

- Use the smallest syringe size possible (e.g., if infusing 2.3 mL of fluid, use a 3 mL syringe)
- Use compatible components which have the smallest internal volume or “dead space”
- Ensure the device is as close to level of the patient’s heart as possible. Patient’s heart should be in line with [CHANNEL SELECT] key
- Use the [PRIME SET WITH SYRINGE] channel option on the Alaris™ Syringe module to speed up the engagement of the device’s mechanical components and decrease the syringe’s internal friction
- If utilizing a pre-run infusion practice (to allow for medication equilibration prior to connection to the patient), ensure the distal end of the administration set is level with or higher than the device
- Avoid use of manifolds with ports containing high pressure valves. These valves require at least 50-200 mmHg pressure to open and allow fluid flow. These high-pressure valves may cause a significant delay in therapy followed by a sudden bolus once the valve is opened, particularly at low infusion rates
- **Note: These recommendations are especially important when infusing high-risk or life- sustaining medications at low infusion rates (for example, <5mL/h and especially at flow rates <0.5mL/h)**

Alaris™ PCA module

Demonstrate PCA Module Setup

- The patient’s heart level should be in line with [CHANNEL SELECT] key.
- System Start Up and Security key lock feature.
 - Use of the security key or security code when installing a new syringe or changing the syringe.
- Priming the set using the “Prime Set with Syringe” feature
 - Prime set prior to attaching to patient.
 - The tubing should be clamped to prevent inadvertent or uncontrolled flow with a primed administration set.
- Loading the syringe into the Alaris™ PCA module
 - Hold the installed syringe plunger to prevent accidental push on the plunger when lowering the drive head.
 - Clamp off fluid flow to the patient before loading and unloading a syringe.
 - Check that the installed syringe matches the manufacturer and size displayed on the pump.

Demonstrate Programming the Alaris™ PCA module

- Programing the following
 - PCA dose + Continuous dose infusion

BD Alaris IV Infusion System Skills Checklist #DAHS-NSCBD18-ALARIS

Page 3 of 4

| | | |
|------------------|---|-----------------------|
| Name: | | Employee ID #: |
| Unit: | | Title: |
| Due Date: | New hire: prior to end of unit orientation period: ____ / ____ / ____. | |
| | Current Staff: | |

These skills will be considered complete when all below performance criteria are completed and pages 1, 2, 3 and 4 have been scanned and emailed to: hs-cppn@ucdavis.edu

| Alaris™ PCA module (<i>Continued</i>) | Date Completed (or N/A) | Verifier Initials |
|---|-------------------------|-------------------|
|---|-------------------------|-------------------|

- How to modify PCA parameters during an active PCA infusion (PCA dose, Lockout interval, Continuous dose, Maximum limit).

Demonstrate Accessing Patient History and the Alaris™ PCA module

- How to view and clear patient history.
- Verbalize that patient history data is stored as a rolling 24-hour time period.
- Verbalize what actions will delete the PCA patient history.

Demonstrate Pausing the infusion, Changing the syringe and Restoring the infusion

- Clamping off fluid flow to the patient before loading and unloading a syringe.
- Pause the infusion, change the current syringe, and then use the [RESTORE] key to restore the previous programming parameters.
- Verbalize that the [RESTORE] key should only be used if the Drug, Therapy, Concentration and Dosing Units remain the same.

Demonstrate Understanding of the Alaris™ PCA Pause Protocol

The Alaris™ PCA module will pause when hospital-established parameters on the Alaris™ etCO₂ module are met.

Demonstrate Understanding of the near end of infusion (NEOI) option.

- Near end of infusion (NEOI) option allows an alert to sound at a hospital-established remaining syringe volume before the infusion is complete (Empty Syringe alert).
- An audio prompt will sound at NEOI, which requires being silenced just once, and will not re-occur following the initial silencing until the empty syringe alert sounds.

I am not responsible for the PCA module.

| Alaris™ EtCO ₂ module | Date Completed (or N/A) | Verifier Initials |
|----------------------------------|-------------------------|-------------------|
|----------------------------------|-------------------------|-------------------|

EtCO₂ Monitoring

- Understanding of the basic parameters monitored using the Alaris™ etCO₂ module, including basic Ventilation vs. Oxygenation and a normal etCO₂ waveform.
- Locating the Gas Inlet on the Alaris™ etCO₂ module and attach the disposable.
Using the directions for use insert as a reference before attaching the disposable to the patient.

Alarms and Limits

- How to view etCO₂ alarm limits, RR alarm limits, and etCO₂ waveform from the main display.
- How change etCO₂, RR, and No Breath limits.

BD Alaris IV Infusion System Skills Checklist #DAHS-NSCBD18-ALARIS

Page 4 of 4

| | | |
|------------------|---|-----------------------|
| Name: | | Employee ID #: |
| Unit: | | Title: |
| Due Date: | New hire: prior to end of unit orientation period: ____ / ____ / ____. | |
| | Current Staff: | |

These skills will be considered complete when all below performance criteria are completed and pages 1, 2, 3 and 4 have been scanned and emailed to: hs-cppn@ucdavis.edu

| Alaris™ EtCO ₂ module (Continued) | Date Completed (or N/A) | Verifier Initials |
|--|-------------------------|-------------------|
|--|-------------------------|-------------------|

Demonstrate Pre-silencing Alarm

- Understands that this mode will only pre-silence the monitoring alarm for 2 minutes and will not silence infusion alarms.

Demonstrate Viewing EtCO₂ Trend Data

- Understand how to view the trend data.
 - How to tell which value has triggered an alarm (bell icon).
 - If there is no data for time period displayed, dashes (---) will be displayed.
- Current patient data will not be displayed while Trend Data feature is being viewed

Demonstrate Understanding of Alarms/Alerts/Troubleshooting

- Verbalize meaning and response to:
 - Auto zero in progress Alarm
 - Disposable Disconnected Alarm
 - Clearing Disposable Alarm
 - Disconnect Occluded Disposable Alarm
- Verbalize possible causes and possible actions to:
 - Low etCO₂ Alarm
 - High etCO₂ Alarm
 - High FiCO₂ Alarm
 - No Breath Detected Alarm

I am not responsible for the etCO₂ module.

PRECEPTOR SIGNATURE

Signature and Printed Name of Preceptor or other verified personnel who have initialed on this form:

| | | |
|-----------------|--------------------|-------------------|
| Initial: | Print Name: | Signature: |
|-----------------|--------------------|-------------------|

PRECEPTEE STATEMENT AND SIGNATURE:

I have read and understand the appropriate UC Davis Health Policies/Procedures and/or equipment operations manual, I have demonstrated the ability to perform the verified skills as noted, and I have the knowledge of the resources available to answer questions.

| Printed Name | Signature | Date |
|--------------|-----------|------|
|--------------|-----------|------|