

<b>Zoll R Series ALS For NICU and Peds Acute Care Checklist # DAHS-NSCNEORSALS19</b>	
Page 1 of 2	
<b>Name:</b>	<b>Employee ID #:</b>
<b>Unit:</b>	<b>Title:</b>
PERFORMANCE CRITERIA - Unless otherwise specified all skills will be demonstrated in accordance with the appropriate UC Davis Health Policy and Procedure.	
These skills will be considered complete when all below performance criteria are completed and pages 1 and 2 are scanned and emailed to <a href="mailto:hs-cppn@ucdavis.edu">hs-cppn@ucdavis.edu</a>	

<b>References:</b>	<b>COMPLETED</b>
1. <a href="#">UC Davis Health Policy 6006: Responding to Medical Emergency Situations (Including Code Blue)</a>	
2. <a href="#">UC Davis Health Policy 6005: Automated External Defibrillator (AED-Zoll)</a>	
Completed the assigned ZOLL R Series ALS Defibrillator tutorials in UC Learning	<input type="checkbox"/>
<b>TEST MODE</b>	
Successfully demonstrates manual 30 Joule defibrillator test, with an understanding that the ZOLL R Series performs self-tests every 24 hours and must remain plugged into the OneStep pads to self-test	<input type="checkbox"/>
Weekly test to be performed by RN	<input type="checkbox"/>
Can check and change paper	<input type="checkbox"/>
<b>MANUAL DEFIBRILLATION</b>	
Locates multifunction cable	<input type="checkbox"/>
Confirms shockable rhythm	<input type="checkbox"/>
Selects defibrillator mode (red)	<input type="checkbox"/>
Presses Energy Select or Charge button	<input type="checkbox"/>
Tells everyone to stand clear	<input type="checkbox"/>
Delivers shock at desired energy level	<input type="checkbox"/>
Defines and adjusts energy levels for Pediatrics (2-4J/kg, up to 10J/kg)	<input type="checkbox"/>
<b>SYNCHRONIZED CARDIOVERSION</b>	
Puts device into SYNC mode.	<input type="checkbox"/>
Selects desired energy (0.5-1 J/kg, up to 2J/kg)	<input type="checkbox"/>
Presses charge button	<input type="checkbox"/>
Tells everyone to stand clear	<input type="checkbox"/>
Delivers synchronized shock	<input type="checkbox"/>
States and demonstrates that SYNC must be activated for each and every synchronous cardioversion	<input type="checkbox"/>

**Zoll R Series ALS For NICU and Peds Acute Care Checklist # DAHS-NSCNEORSALS19**

Page 2 of 2

<b>Name:</b>	<b>PPS#:</b>
<b>Unit:</b>	<b>Title:</b>
PERFORMANCE CRITERIA - Unless otherwise specified all skills will be demonstrated in accordance with the appropriate UC Davis Health Policy and Procedure.	
These skills will be considered complete when all below performance criteria are completed and checklist is scanned and emailed to <a href="mailto:hs-cppn@ucdavis.edu">hs-cppn@ucdavis.edu</a>	

<b>PACING</b>	
Turns up pacing output (mA) until capture is achieved – identifies capture	<input type="checkbox"/>
Adjusts pace rate	<input type="checkbox"/>
Understands pausing for visualizing patients underlying rhythm	<input type="checkbox"/>
<b>PADS</b>	
Connects OneStep™ pads to OneStep cable (or other pads/paddles if applicable)	<input type="checkbox"/>
Opens OneStep packaging correctly	<input type="checkbox"/>
Demonstrates proper pad placement for defibrillation, pacing, and cardioversion	<input type="checkbox"/>
Identifies CPR Sensor and explains its purpose	<input type="checkbox"/>
<b>SUPERUSER/TRAINER</b>	
Can access data from the code (Print Chart, Print Log, or Transfer Data)	<input type="checkbox"/>
User demonstrates sufficient understanding of device to train other users in its use	<input type="checkbox"/>

<b>Preceptor Signature: Signature and Printed Name of Preceptor or other verified personnel who have initialed on this form:</b>		
Initial:	Print Name:	Signature:

**PRECEPTEE STATEMENT AND SIGNATURE:**

I have read and understand the appropriate UC Davis Health Policies/Procedures and/or equipment operations manual, I have demonstrated the ability to perform the verified skills as noted, and I have the knowledge of the resources available to answer questions.

Print Name	Date	Signature