

<b>ZOLL X Series Skills Checklist #DAHS-NSCX17</b>	
Page 1 of 3	
Name:	Employee ID #:
Unit:	Title:
PERFORMANCE CRITERIA - Unless otherwise specified all skills will be demonstrated in accordance with the appropriate UC Davis Health Policy and Procedure.	
These skills will be considered complete when pages 1, 2 and 3 are completed and scanned and emailed to <a href="mailto:hs-cppn@ucdavis.edu">hs-cppn@ucdavis.edu</a>	

References:	Date Completed (or N/A)	Verifier Initials
<a href="#">UC Davis Health Policy 6005</a> : Automated External Defibrillator (AED-Zoll)		
Completed the assigned ZOLL X Series Defibrillator tutorials in UC Learning		
<b>TEST MODE</b>		
Successfully demonstrates 30 Joule defibrillator test		
Can check and change paper		
<b>MONITOR MODE</b>		
Applies 3-lead or 12-lead ECG		
Changes lead (ECG Tracing) views		
Locates Recorder key and prints a strip		
Access HR menu and demonstrate how to change settings		
Locate NIBP soft key and activate manual BP measurement		
Access NIBP menu and verbalize options		
Demonstrate how to change NIBP alarm settings		
Change NIBP mode from Manual to Automatic		
Change Automatic mode intervals		
Access SpO2 menu and verbalize options		
Access Temp menu and verbalize options		
Locate CO2 soft key and press to activate		
Access CO2 menu and verbalize options		
Access IBP menu and verbalize options		
<b>MANUAL DEFIBRILLATION</b>		
Locates multifunction cable		
Confirms shockable rhythm		
Selects defibrillator mode (red)		
Presses Energy Select or Charge button		
Tells everyone to stand clear		
Delivers shock at desired energy level		
Defines and adjusts energy levels for Adults (120,150, 200J) and Pediatrics (2-4J/kg)		

**ZOLL X Series Skills Checklist #DAHS-NSCXS17**

Page 2 of 3

<b>Name:</b>	<b>PPS#:</b>
<b>Unit:</b>	<b>Title:</b>
PERFORMANCE CRITERIA - Unless otherwise specified all skills will be demonstrated in accordance with the appropriate UC Davis Health Policy and Procedure.	
<b>These skills will be considered complete when pages 1, 2 and 3 are completed and scanned and emailed to <a href="mailto:hs-cppn@ucdavis.edu">hs-cppn@ucdavis.edu</a></b>	

<b>CPR FEEDBACK</b>	<b>Date Completed (or N/A)</b>	<b>Verifier Initials</b>
Demonstrates steps to fill CPR Index™ – understands proper rate/depth		
Shows that if rate is too slow, metronome beeps and Rate prompt appears		
Speeds up to silence metronome and allow the Rate prompt to disappear		
Shows that if depth is too shallow, the Depth prompt appears on the screen		
Pushes hard to allow Depth prompt to disappear		
Demonstrates understanding of See-Thru CPR® filtered ECG		
<b>SYNCHRONIZED CARDIOVERSION</b>		
Puts device into SYNC mode		
Selects desired energy. (Pediatric stat with 0.5-1.0 J/kg)		
Presses charge button		
Tells everyone to stand clear		
Delivers synchronized shock		
States and demonstrates that SYNC must be activated for each and every synchronous cardioversion		
<b>PACING</b>		
Turns up pacing output (mA) until capture is achieved – identifies capture		
Adjusts pace rate		
Understands pausing for visualizing patients underlying rhythm		
<b>PADS</b>		
Connects OneStep™ pads to OneStep cable (or other pads/paddles if applicable)		
Opens OneStep packaging correctly		
Demonstrates proper pad placement for defibrillation, pacing, and cardioversion		
Identifies CPR Sensor and explains its purpose		
<b>INTERNAL PADDLES</b>		
Understands how to connect internal paddles to OneStep™ cable		
Selects defibrillator mode (red)		
Understands 10J default energy level with range of 1 to 50 Joules		

**ZOLL X Series Skills Checklist #DAHS-NSCXS17**

Page 3 of 3

<b>Name:</b>	<b>PPS#:</b>
<b>Unit:</b>	<b>Title:</b>
PERFORMANCE CRITERIA - Unless otherwise specified all skills will be demonstrated in accordance with the appropriate UC Davis Health Policy and Procedure.	
<b>These skills will be considered complete when pages 1, 2 and 3 are completed and scanned and emailed to <a href="mailto:hs-cppn@ucdavis.edu">hs-cppn@ucdavis.edu</a></b>	

SUPERUSER/TRAINER	Date Completed (or N/A)	Verifier Initials
Demonstrate how to use additional options (Mentor mode, Set the clock, etc.)		
Understands how to change parameter settings (NIBP, EtCO2, SpO2)		
Understands purpose of Code Marker		
Can access data from the code (Print Chart, Print Log, or Transfer Data)		
User demonstrates sufficient understanding of device to train other users in its use		

<b>Preceptor Signature: Signature and Printed Name of Preceptor or other verified personnel who have initialed on this form:</b>		
Initial:	Print Name:	Signature:

**PRECEPTEE STATEMENT AND SIGNATURE:**

I have read and understand the appropriate UC Davis Health Policies/Procedures and/or equipment operations manual, I have demonstrated the ability to perform the verified skills as noted, and I have the knowledge of the resources available to answer questions.

Print Name	Date	Signature