

Hypotension Prediction Index: Correlations between Invasive and Non-invasive Pressure Inputs

Kamron Sarhadi, BS; Jorge Avila, MD; Neal Fleming, MD, PhD

Background


- Continuous BP monitoring is essential to intraoperative care, as hypotensive events can significantly increase the risk of AKI, MI, and mortality post-op^{1,2}.
- The Hypotension Prediction Index (HPI) is a novel algorithm derived from machine learning that gives anesthesiologists the ability to predictive hypotensive events.
- The HPI derived from intra-arterial catheter monitoring has been shown to predict hypotensive events with sensitivity and specificity >80%³.
- The utility and accuracy of the HPI when derived from non-invasive monitoring techniques, such the ClearSight finger cuff, have yet to be examined.
- This study seeks to compare the intra-arterial catheter-derived HPI vs the ClearSight finger cuff-derived HPI (as well as other hemodynamic variables), to see if it is viable tool for anesthesiologists to use when invasive monitoring is not indicated.

Methods

- Recorded patients' hemodynamics concurrently with both invasive (intra-arterial) and non-invasive (ClearSight) monitoring.
- Each monitoring system was connected to a HemoSphere monitor with the HPI software.
- Data collected from the ClearSight system was compared to corresponding intra-arterial waveform data using Pearson correlation analysis, Bland-Altman analysis, and analysis of concordance.

Results


Procedures by Specialty




Demographics

NUMBER OF PATIENTS	AGE	HEIGHT	WEIGHT	SEX
42	66.7 ± 17.2	170.9 ± 10.7	87 ± 22.5	59.5% M, 40.5% F


Raw Data (HPI)




Correlation



r	0.78
95% CI	0.77 to 0.78
r²	0.60




r	0.77
95% CI	0.76 to 0.78
r²	0.59




r	0.72
95% CI	0.71 to 0.73
r²	0.52


Bland-Altman



Bias	+9.8
95% LOA	± 30




Bias	+3.4
95% LOA	± 19




Bias	-11
95% LOA	± 46


Concordance



Percent Concordance	77%
----------------------------	------------



Percent Concordance	78%
----------------------------	------------




Percent Concordance	72%
----------------------------	------------

Discussion

- Correlation between HPI values is strong with r value > 0.7.
- Bland-Altman: Differences between ClearSight vs. Intra-arterial HPI becomes less pronounced at the extremes and more pronounced at values between 35-80.
- Strong correlation & percent concordance between ClearSight and Intra-arterial MAP and Systolic Pressure.
- Overall, the ClearSight finger cuff presents as a promising monitoring system that compares well to the gold standard for hemodynamic monitoring with some drawbacks, including occasionally lacking reliability and producing significant noise.

Next Steps

- Expanding analysis to include full range of 11 hemodynamic variables including CO, SV, SVV, Eadyn, and PPV.
- Error grid analysis



References

1. Salmasi V, Maheshwari K, Yang D, et al. *Anesthesiology*. 2017;126:47-65.
2. Walsh M, Devereaux PJ, Garg AX, et al. *Anesthesiology*. 2013;119:507-515.
3. James Davies S, Tilma Vistisen S, Jian Z, et al. *Perioperative Medicine*. 2020;130:352-359.

Acknowledgements

Thank you to Edwards Lifesciences for their help with data extraction and error grid plot analysis.

Kamron Sarhadi, BS, UC Davis SOM
Jorge Avila, MD, UC Davis SOM
Neal Fleming, MD, PhD, UC Davis SOM
 Department of Anesthesiology and Pain Medicine

Download poster & Abstract:



UCDAVIS
SCHOOL OF MEDICINE